

RAYSTATION v2025

IHE-RO Integration Statement

Safety notices

Warning and caution notices in the user documentation inform about the safe use of the product and must be followed.



WARNING!

A warning notice informs about a risk for bodily harm or death. In most cases, the risk is related to mistreatment of the patient.



CAUTION!

A caution notice informs about a risk for damage to equipment, software or data.

Note: A note provides additional helpful information, tips or reminders.

Copyright

This document contains proprietary information that is protected by copyright. No part of this document may be photocopied, reproduced or translated to another language without prior written consent of RaySearch Laboratories AB (publ).

All Rights Reserved. © 2025, RaySearch Laboratories AB (publ).

Printed material

Hard copies of Instructions for Use and Release Notes related documents are available upon request.

Trademarks

RayAdaptive, RayAnalytics, RayBiology, RayCare, RayCloud, RayCommand, RayData, RayIntelligence, RayMachine, RayOptimizer, RayPACS, RayPlan, RaySearch, RaySearch Laboratories, RayStation, RayStore, RayTreat, RayWorld and the RaySearch Laboratories logotype are trademarks of RaySearch Laboratories AB (publ)*.

Third-party trademarks as used herein are the property of their respective owners, which are not affiliated with RaySearch Laboratories AB (publ).

RaySearch Laboratories AB (publ) including its subsidiaries is hereafter referred to as RaySearch.

* Subject to registration in some markets.

INTEGRATION STATEMENT

1.1 INTRODUCTION

This IHE Integration Statement is prepared and published to describe the conformance of RayStation with the IHE-RO Technical Framework. It identifies specific IHE-RO capabilities that are supported in terms of actors and integration profiles described in the technical frameworks of each domain.

Users familiar with these concepts can use Integration Statements to determine the integration of RayStation with complementary systems and what clinical and operational benefits such integration might provide. This Integration Statements is intended to be used in conjunction with the RayStation DICOM Conformance Statement.

1.2 PRODUCT VERSION

This Integration Statement is valid for RayStation v2025.

1.3 ABOUT THE TESTING PROCESS

IHE provides a process for vendors to test their implementations of IHE actors and integration profiles. The IHE testing process, culminating in a multi-party interactive testing event called the Connect-a-thon, provides vendors with valuable feedback and provides a baseline indication of the conformance of their implementations. The process is not intended to independently evaluate, or ensure, product compliance. In publishing the results of the Connect-a-thon and facilitating access to vendors' IHE Integration Statements, IHE and its sponsoring organizations are in no way attesting to the accuracy or validity of any vendor's IHE Integration Statements or any other claims by vendors regarding their products.

RaySearch Laboratories has done further validation beyond the IHE-RO Connect-a-thon to ensure that RayStation conforms to applicable standards. This Integration Statement shall not be used a guarantee that RayStation will work in any environment or with any external system.

1.4 REFERENCES

Reference	Resource
DICOM	https://www.dicomstandard.org/
IHE	https://www.ihe.net/
IHE-RO	https://www.ihe.net/ihe_domains/radiation_oncology/

1.5 DEFINITIONS

Term	Meaning
CE	CE marking is a certification mark that indicates conformity with health, safety, and environmental protection standards for products sold within the European Economic Area (EEA)
DICOM	Digital Imaging and Communications in Medicine
EEC	European Economic Community
IHE	Integrating the Healthcare Enterprise
IHE-RO	Integrating the Healthcare Enterprise Radiation Oncology
BRT0-II	Basic RT Object Interoperability-II (IHE-RO Profile)
DRRO	Deformable Registration for Radiation Oncology (IHE-RO Profile)
MMRO-III	Multi-Modality Registration in Radiation Oncology-III (IHE-RO Profile)
TPPC	Treatment Planning – Plan Content (IHE-RO Profile)

1.6 INTEGRATION STATEMENT

IHE-Integration Statement	Date	5 December 2025
Vendor	Product Name	Version
RaySearch Laboratories	RayStation	v2025
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
BRT0-II	Contourer	Resample CT, Feet First, Structure Set Retrieval
	Dose Displayer	-
MMR0-III	Registrator	-
	Registered Contourer	-
	Registered Display	-
	Registered Dose Display	-
DRR0	Deformable Registrator	-
	Deformable Displayer	-
	Image Deformer	-
	Dose Deformer	-
	Contour Deformer	-
TPPC	Basic Static Beam Producer	Bolus and Block Beam Modifier
	Basic Static MLC Beam Producer	Bolus and Block Beam Modifier
	Arc Beam Producer	Bolus Beam Modifier
	MLC Fixed Aperture Arc Beam Producer	Bolus Beam Modifier
	MLC Variable Aperture Arc Beam Producer	Bolus Beam Modifier
	Hard Wedge Beam Producer	Bolus Beam Modifier
	Virtual Wedge Beam Producer	Bolus Beam Modifier
	Motorized Wedge Beam Producer	Bolus Beam Modifier
	Static Electron Beam Producer	Bolus and Block Beam Modifier
	Step & Shoot Beam Producer	Bolus Beam Modifier
	Sliding Window Beam Producer	Bolus Beam Modifier
	IMAT/VMAT Beam Producer	Bolus Beam Modifier
	Photon Applicator Beam Producer	Bolus Beam Modifier
	Photon Applicator Arc Beam Producer	Bolus Beam Modifier
	Basic Static Beam Consumer	Bolus Beam Modifier
	Basic Static MLC Beam Consumer	Bolus Beam Modifier
	Arc Beam Consumer	Bolus Beam Modifier

	MLC Fixed Aperture Arc Beam Consumer	Bolus Beam Modifier
	MLC Variable Aperture Arc Beam Consumer	Bolus Beam Modifier
	Hard Wedge Beam Consumer	Bolus Beam Modifier
	Virtual Wedge Beam Consumer	Bolus Beam Modifier
	Motorized Wedge Beam Consumer	Bolus Beam Modifier
	Static Electron Beam Consumer	Bolus and Block Beam Modifier
	Step & Shoot Beam Consumer	Bolus Beam Modifier
	Sliding Window Beam Consumer	Bolus Beam Modifier
	IMAT/VMAT Beam Consumer	Bolus Beam Modifier
	Photon Applicator Beam Consumer	Bolus Beam Modifier
	Photon Application Arc Beam Consumer	Bolus Beam Modifier
Links to Standards Conformance Statements for the Implementation		
DICOM IHE	https://www.raysearchlabs.com/product-configurations/	
Links to general information on IHE		
In North America: www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j



CONTACT INFORMATION



RaySearch Laboratories AB (publ)
Eugeniavägen 18C
SE-113 68 Stockholm
Sweden

Contact details head office

P.O. Box 45169
SE-104 30 Stockholm, Sweden
Phone: +46 8 510 530 00
Fax: +46 8 510 530 30
info@raysearchlabs.com
www.raysearchlabs.com

RaySearch Americas

Phone: +1 347 477 1935

RaySearch China

Phone: +86 137 0111 5932

RaySearch India

E-mail: manish.jaiswal@raysearchlabs.com

RaySearch Singapore

Phone: +65 8181 6082

RaySearch Australia

Phone: +61 411 534 316

RaySearch France

Phone: +33 (0)1 76 53 72 02

RaySearch Japan

Phone: +81 (0)3 44 05 69 02

RaySearch UK

Phone: +44 (0)2039 076791

RaySearch Belgium

Phone: +32 475 36 80 07

RaySearch Germany

Phone: +49 (0)172 7660837

RaySearch Korea

Phone: +82 01 9492 6432